

BIOKINEKT




Neurocor
trust motion

BALANCE SYSTEM



new technologies – new solutions

Biokinekt with Balance system «Romberg» and biomechanical sensors «Bisens» allows to conduct in-depth diagnostics and highly effective rehabilitation using biofeedback.



BALANCE SYSTEM


METHODS:



Classical stabilometry on the support reaction



Controlling myograms of various muscle groups



3D Stabilometry and analysis of different human body segments tremor.

Fascinating gaming environments with the possibility of «immersion» to the virtual reality provides increased patient motivation, the treatment process becomes pleasant and interesting. BLOKINEKT complex allows you to use the following physiological parameters: support reaction (on the coordinates of the Common Central Pressure (CCP) on the plane), movements in segments, rotations in joints and segments, tremor and electrical muscle activity

BALANCE SYSTEM



VR biofeedback programs with a large-screen display or VR helmet are implemented.



Body weight support system is used for the early stage of rehabilitation, that allows to rehabilitate people with disabilities.



Neurocor LLC
develops
and manufactures
equipment
for medical
rehabilitation
and diagnostics
uses in health
care and sports.
We are open
for cooperation
and always glad
to welcome partners
and associates.



Contact
our
manager



Neurocor
trust motion

 **Made in Russian Federation**

Neurocor LLC
www.neurocor.uk
info@neurocor.ru
+7 (495) 740-35-34

109428, Russian Federation, Moscow, 24/1 Ryazanskiy prospekt